


- These alphabetical codes are taken from the diagrams on the next page.

	A			C			F			G			E		
Size	Slim	Med	Wide	Slim	Med	Wide	Slim	Med	Wide	Slim	Med	Wide	Slim	Med	Wide
20	15.4		18.5	21.2		25.6	5.7		6.6	4.0		4.9	13.3	13.3	13.3
21	15.8		18.9	22.0		26.3	5.8		6.7	4.1		5.1	14.0	14.0	14.0
22	16.2		19.3	22.8		27.0	6.0		6.8	4.2		5.2	14.7	14.7	14.7
23	16.6		19.7	23.6		27.7	6.1		7.1	4.3		5.4	15.3	15.3	15.3
24	17.0	19.5	20.1	24.4	27.8	28.4	6.2	6.8	7.2	4.4	4.9	5.5	16.0	16.1	16.0
25	17.4	19.9	20.5	25.2	28.5	29.1	6.3	6.9	7.3	4.4	5.0	5.6	16.7	16.8	16.7
26	18.2	20.2	20.9	26.5	29.2	29.8	6.5	7.0	7.5	4.4	5.1	5.7	17.3	17.4	17.3
27	18.5	20.9	21.3	27.2	30.0	30.5	6.6	7.1	7.6	4.6	5.2	5.9	18.0	18.0	18.0
28	18.9	21.1	21.8	28.0	30.7	31.4	6.7	7.4	7.9	4.8	5.3	6.0	18.7	18.7	18.7
29	19.3	21.6	22.2	28.8	31.4	32.1	6.8	7.3	8.0	4.9	5.4	6.1	19.3	19.4	19.3
30	19.8	22.0	22.7	29.5	32.2	32.8	7.0	7.5	8.1	4.9	5.5	6.1	20.0	20.1	20.0
31	20.3	22.4	23.2	31.0	32.9	33.5	7.2	7.7	8.6	4.9	5.6	6.2	20.7	20.7	20.7
32	20.7	22.8	23.6	31.8	33.6	34.2	7.4	7.8	8.7	5.0	5.7	6.2	21.3	21.4	21.3
33	21.1	23.1	23.8	32.6	34.3	34.9	7.5	7.9	9.0	5.1	5.8	6.4	22.0	22.0	22.0
34	21.6	23.7	24.4	33.3	35.4	35.6	7.6	8.1	9.1	5.1	5.9	6.4	22.7	22.7	22.7
35	22.0	24.2	24.8	34.0	36.2	36.3	7.8	8.2	9.3	5.3	6.0	6.5	23.3	23.3	23.3
36	22.4	24.6	25.2	34.8	37.0	37.0	7.9	8.4	9.3	5.4	6.1	6.5	24.0	24.0	24.0
37	22.7	25.2	25.6	35.6	37.5	37.7	8.0	8.5	9.4	5.7	6.1	6.6	24.7	24.7	24.7
38	22.4	25.5	26.0	35.6	38.2	38.4	8.2	8.7	9.6	5.8	6.2	6.6	25.3	25.3	25.3
39	22.8	25.8	26.4	37.4	39.0	39.1	8.4	8.8	9.8	5.9	6.3	6.7	26.0	26.0	26.0
40	23.2	26.2	26.8	38.4	39.8	39.8	8.5	9.0	9.9	6.0	6.4	6.8	26.7	26.7	26.7
41	23.6	26.7	27.2	39.2	40.4	40.5	8.6	9.1	9.9	6.2	6.5	6.8	27.3	27.3	27.3
42	24.0	27.1	27.6	39.8	40.8	41.2	8.9	9.2	10.1	6.2	6.6	7.0	28.0	28.0	28.0

NOTE!

The table contains the minimum measurements of the orthosis and is used solely to assess whether the orthosis will fit in the shoe. The measurement of the heel width and forefoot width is final. The margin of adjustment occurs only in the circumference measurement. Due to the adjustment of the lace and the softness of the leather, the circumferences of the orthosis may be slightly larger than the given minimum. Standard footwear may not be perfectly fitted to the orthosis. The straps in the footwear must be adjusted and shortened by yourself. Any non-standard orthoses require individual implementation. To prevent the footwear from being too tight, add about 1.5 cm to the measurements of the insole length.

- A. Circumference of forefoot/ at the MTP joints
- B. Circumference of midfoot / at highest point of the arch
- C. Circumference around the heel and ankle / back of heel to front of ankle
- D. Circumference above the ankle / above the malleoli



- E. Length of the foot
- F. Width of the forefoot at the MTP joints
- G. Width of the heel

