



DATE MEASUREMENTS TAKEN

Participant Details

Name:	<input type="text"/>		
Date of Birth:	<input type="text"/>	<input type="checkbox"/> Male	<input type="checkbox"/> Female
Diagnosis/ Disability:	<input type="text"/>		
NDIS Number:	<input type="text"/>	<input type="checkbox"/> Self Managed	<input type="checkbox"/> Plan Managed
Plan Manager invoice email:	<input type="text"/>		

Contact Details

Parent Name:	<input type="text"/>
Phone Number:	<input type="text"/>
Email:	<input type="text"/>
Address:	<input type="text"/>

Details of prescribing Allied Health Professional

Name:	<input type="text"/>
Allied Health Profession :	<input type="text"/>
Organisation :	<input type="text"/>
Email :	<input type="text"/>
Delivery Address (if different to address in contact details):	<input type="text"/>

Made-to-Measure Hensinger Collar

This form is used to place an order for a made-to-measure Hensinger Collar. Complete pages 1 and 2 of the form and return via email. Please be aware that made-to-measure Hensinger collars cannot be returned or refunded.

Signature (if submitted via email this is taken as an electronic agreement to this declaration)

By submitting this form, I agree that I, as a registered health professional, have assessed the client and prescribed the Hensinger collar with the specifications in this order with due diligence and care and that the collar is non-returnable and non-refundable. I agree for the measurements and any photographs attached to this document to be shared with Danmar Products in the USA for manufacturing purposes.

Name: _____ Signed: _____



Made-To-Measure Hensinger Collar Measurements


For accurate measurements for a made-to-measure Hensinger collar, consider the following details:

- Ensure head and neck is supported in a neutral alignment and breathing, swallowing and vision functions are maintained or improved.
- When using a chest support, subtract its thickness from the "B" measurements. Straps or trunk supports of 3mm or less should not affect the fit, while anything thicker than 3mm should be taken into consideration.
- Take into account the location of gastrostomy sites when measuring the length of the collar ears.

A - Circumference

B - Height

C - Width



Hensinger Type:

☐

Low back

☐

High back

	Made-To-Measure Hensinger Collar Dimensions	Measurements in cm	Office Use Only: Inches
A	Circumference (measure neck circumference, round to nearest 1cm)		
B	Height of the collar (measure vertically from top of shoulder to bottom of earlobe, then deduct 12mm)		
C	Width of the collar		
D	Length of collar ears (measure from chin to the upper 1/4 - 1/3 of sternum with head and neck in a neutral position)		

Any specific requirements or relevant measurements (attach photographs of the head and neck and the specific issue that requires addressing):

Colour of Lycra covering (one standard cover is supplied with the collar):

Number and colour of additional Lycra coverings required (purchased in a pack of 2 covers):

