

Straight and Reverse Last Footwear



Straight Last

The footwear is shaped neutrally, with neither forefoot abduction or adduction. The heel counter is high and firm to maintain rearfoot alignment, whilst the medial stiffener is elongated to assist with maintaining a corrected neutral forefoot position.



Reverse Last

The footwear is shaped with forefoot abduction. The heel counter is high and firm to maintain rearfoot alignment, whilst the medial stiffener is elongated to assist with maintaining an abducted forefoot position.

Style options available from AurelkaOrto in both straight and reverse last are shown below. All styles have several colour options. Please contact Well & Able for colour options once the last and style have been selected.



AurelkaOrto 1001 Boot



AurelkaOrto 1011 Boot



AurelkaOrto 1003 Sandal



AurelkaOrto 1016 Closed-toe Sandal