

my generation



by **schein**



Orthotic shoes with

BOA®

fastening system

Shoe care

Orthotic shoes

Space for therapy

Many people find wearing a leg orthosis too restricting for their freedom of movement. To combat this, our orthotic shoes are designed to ensure an optimum fit for the foot, orthosis and therapeutic shoe. The sturdy heel cap combined with the flexible raised collar guarantees a good hold and a secure stance.


However, an ideal fit is not the only thing that matters for our shoes. For ease of everyday use, it is also vitally important that our orthotic shoes are easy to put on and take off. The Boa® fastening system makes fastening and unfastening with one hand possible.



However, when designing our orthotic shoes, we do not confine ourselves to addressing the medical issues: fashion considerations also play an important role.

That is why our therapeutic shoes strongly resembles normal shoes in design. To ensure this, we keep track of the latest fashions, popular cuts and colours and integrate these into our therapeutic shoes to help people feel good in them.

If further adjustments are needed for orthotic purposes, our orthopaedic shoe technicians can make these on an individual basis. Specific orthopaedic adaptations and fittings can be added to all our orthotic shoes, ensuring that every foot receives optimum care.



Available in the widths TN5, TN8, TN9 and TN12. Where each foot requires different care, shoes of differing widths and sizes can be combined without difficulty.

Illustration: diagrammatic depiction of insole shape and internal volume of the last.

Orthotic shoes

Our orthotic shoes:
a combination of proven
medical functionality and a
fresh appearance

The right footwear is an important addition to an individually produced orthosis. Shoes worn by people wearing an orthosis not only have to satisfy certain key requirements, but must also make walking and running easier.

Due to their material quality and shape, conventionally manufactured shoes are not up to meeting the requirements of orthoses. That is why we recommend our therapeutic shoes.

Flexible outsole

Ever thinner and more flexible materials are used in the orthopaedic handcraft and its treatments.

In response, we have introduced an even more flexible version of the sole for our orthotic shoes.

All shoes designated in the tables as 'flexible' are extremely elastic. None of this has any negative impact on the therapeutic benefits of the orthosis.



Taken care of with a single turn



Orthotic shoes with **Boa®** fastening system

Comfortable to put on
thanks to the flexible
hinged back

Even tightening of
the nylon cord
ensures ideal hold

Sturdy
heel cap

Sturdy toe cap
to protect the
toes



Practical fastening and
width adjustment with
Boa® fastening system – either
on the side at the front or on the back

Our therapeutic shoes are specially designed for people who have to wear an orthosis. The combination of specialist orthopaedic knowledge, modern leathers and high-quality materials makes for long-lasting companions.

The Boa® fastening system supports the philosophy of our orthotic shoes thanks to its robust workmanship:

- Very good hold
- Fine adjustment
- One-handed operation
- Low weight
- Available in two widths
- Available as a front or rear fastening

Boa® fastening system

We offer fastening with a twist

We have fitted all the orthotic shoe models in this brochure with a Boa® fastening system. Originally developed for use in sport, this quick-lacing system for wearers of orthotics with limited motor functions makes putting on and taking off the shoes considerably easier. This is because the fastening and unfastening is done with a single movement of the adjustment wheel. Turning the wheel adjusts the entire lacing evenly and individually for the patient and fixes it to the millimetre.

The fastening ensures a stable hold and the highest levels of comfort when wearing the shoe.

The secure hold

We offer the Boa® fastening system in two versions:


- With adjustment wheel on the back
- With adjustment wheel on the outer side

Both lacing systems give the foot excellent hold by achieving even tightening across the entire lace surface.


Unlike conventional shoelaces, which can absorb water or dirt, the low-friction synthetic coverings attached to the fine stainless-steel wires remain feather-light.

Even under heavy strain, the cords and wires are a long-lasting concept.

Once it has been fastened, the lacing stays in place – the orthosis and/or foot are securely fixed.



To unfasten:
just pull out the
adjustment wheel
and the fastening is
released.



To fasten:
push down the
adjustment wheel and
then set the fastening
by turning the wheel
clockwise





Flexible
sole

Applies to all models
on this double page

Merkur

362910 / width TN5
372910 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole
Uranus	362900 372900	TN5 TN8	28 - 42	Nubuck	black / grey	Boa® fastening system	textile	flexible
Merkur	362910 372910	TN5 TN8	28 - 35	Nubuck	purple / grey	Boa® fastening system	textile	flexible
Neptun	362920 372920	TN5 TN8	28 - 42	Nappa / Nubuck	blue / grey	Boa® fastening system	textile	flexible

Orthotic shoes with Boa® fastening system



Uranus

362900 / width TN5
372900 / width TN8



Neptun

362920 / width TN5
372920 / width TN8





Moros

362300 / width TN5
372300 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining
Moros	362300 372300	TN5 TN8	30 - 42	Nappa / Mesh fabric	black	Boa® fastening system	textile
Leto	362310 372310	TN5 TN8	24 - 42	Nubuck / Mesh fabric	grey / dark blue	Boa® fastening system	textile
Hera	362320 372320	TN5 TN8	24 - 32	Nubuck / Mesh fabric	pink	Boa® fastening system	textile
Hesti	362330 372330	TN5 TN8	30 - 42	Nappa / Faded denim look	beige	Boa® fastening system	textile

Orthotic shoes with Boa® fastening system



Leto

362310 / width TN5
372310 / width TN8



Hesti

362330 / width TN5
372330 / width TN8



Hera

362320 / width TN5
372320 / width TN8



my design CONFIGURATOR

Now also with



My individual shoe


You can put together your own unique pair of shoes just how you want them with our my design shoe configurator.

From the sole to the inner lining – and now also with the Boa® fastening system – there are numerous design possibilities.



Configuration examples
We reserve the right to make changes to colour and designs

1



Configuration code:
DESBSC01




Configuration code:
DESBSC02



Configuration code:
DESBSC03

2



Configuration code:
DESBSP01



Configuration code:
DESBSP02



Configuration code:
DESBSP03


	Model	Item no.	Width	Size	Fastener
Orthotic shoe	1	362930 372930	TN5 TN8	28 - 42	Boa® fastening system
Orthotic shoe	2	362340 372340	TN5 TN8	24 - 42	Boa® fastening system

It's this simple

1. Select model



2. Select pair/individual and size/width



3. Design to suit your wishes



4. Download order as a PDF



Further orthotic shoes and much more in our catalogue

Available from our
my generation partners
or at
mygeneration-schein.de



my generation

by **schein**



We reserve the right to make minor changes to colours and models.
Illustrated shoes are size 28.

Your specialist dealer:



Schein Orthopädie Service KG
Hildegardstr. 5
42897 Remscheid
Germany
Tel. +49 2191 910-201
Fax +49 2191 910-102
remscheid@schein.de
www.mygeneration-schein.de

carbon neutral
natureOffice.com | DE-221-517441
print production

schein orthopädie
service
Germany since 1879