



Orthotic • Support • Anti-varus
Shoes



Shoe collection

Orthotic • Support • Anti-varus
Shoes



www.mygeneration-schein.de



Foreword



A great feeling, and not just for the foot

For four generations, our company has stood for a conscious and sensitive approach to orthopaedic (shoe) technology, people and mobility. Our more than 30 years of experience in the production of therapeutic shoes and the knowledge passed down from generation to generation guarantee high quality and the best workmanship when it comes to our shoes.

The philosophy underlying every single therapeutic shoe we design is to take the principles and findings of medical science as our benchmark for ensuring that we provide the best possible care for your feet. We are continuously working on the further development of our products and are in constant dialogue with our customers and doctors.

We incorporate fashionable elements, based on the latest trends and colours, into our shoes. This enables us to produce orthopaedic shoes that meet our aim of resembling normal shoes in appearance. More than 3,000 partners across the world carry our my generation products, guaranteeing you international availability.

My generation stands for high-quality and environmentally friendly manufacture combined with a love of detail, a high proportion of craftsmanship and the use of high-quality materials. In this way, we make our contribution towards ensuring the best possible foot care for you. And this is something you can rely on.

A handwritten signature in blue ink, appearing to read "Arnd Schein".

Arnd Schein
Schein Orthopädie Service KG

Table of contents

Sizing	8
Shoe care	10–17
Insole care	18–19
Quality management	20
Materials	21–22
Sole types	24
Sole materials, heels, etc.	26–28
The Schein app	29

my design configurator

my design configurator and models	30–35
---	-------

Orthotic shoes

36–91

Support shoes

92–111

Multi-support shoes

112–113

Anti-varus shoes

114–143





Sizing

The shoe must fit

Our feet have to bear our weight all through our lives. This means our shoes do a lot of important jobs. They have to support us without us noticing them. They must be durable and comfortable. And they must not have a negative impact on growth.

The foot of a small child up to the age of four generally grows at a rate of up to two centimetres a year, which is up to three shoe sizes. From the age of four onwards, annual growth can be up to two shoe sizes. Because of this, shoes should never be bought too small. They also mustn't be too big, because overlarge shoes can harm the foot as much as ones that are too small.

Our therapeutic shoes ensure the right fit

A good heel fit plus versatile adjustment options across the instep ensure a well-fitting shoe. Our shoes offer ample space for the toes in the forefoot region.

Working out the right size

For children and teenagers, the shoe size is made up as follows:

Foot length + push room + room to grow = shoe size
 Allowance (add-on)

We have designed a special foot gauge for our shoes, which takes into account the necessary allowance. On the back of the foot gauge is a table showing the widths of our orthotic and support shoe range.



Item no.: 032640-002



Orthotic shoes

Space for therapy

Many people find wearing a leg orthosis too restricting for their freedom of movement. To combat this, our orthotic shoes are designed to ensure an optimum fit for the foot, orthosis and therapeutic shoe. The sturdy heel cap combined with the flexible raised collar guarantee a good hold and a secure stance.



However, an ideal fit is not the only thing that matters for our shoes. For ease of everyday use, it is also vitally important that our orthotic shoes are easy to put on and take off. Therefore many

of our models feature a slip-in loop which makes it easier to do just that. But that is not all: the special 'Derby' construction allows the whole shoe to be opened up very wide.

If further adjustments are needed for orthotic purposes, our orthopaedic shoe technicians can make these on an individual basis. Specific orthopaedic adaptations and fittings can be added to all our orthotic shoes, ensuring that every foot receives optimum care.



Great all-round
foot care



Available in the widths TN5, TN8, TN9 and TN12. Where each foot requires different care, shoes of differing widths and sizes can be combined.

Illustration: diagrammatic depiction of insole shape and internal volume of the last.

Support shoes

The right footwear provides valuable support


Our support shoes are suitable for a wide range of uses. They can relieve the strain on the foot and ankle joints as well as provide them with support and guidance.

They are proven for the following indications:

- Support for painful foot joints in cases of rheumatoid arthritis or arthritis for any other cause
- Support for joints with weak ligaments after injuries to ankle joint ligaments or in cases of congenital fibular aplasia
- Guidance for joints which are neurologically impaired due to cerebral palsy, spina bifida or poliomyelitis
- Support for slight to moderate foot deformities



Wide range of fittings
for various different indications



Available in the widths TN2, TN4 and TN6. Where each foot requires different care, shoes of differing widths and sizes can be combined.

Illustration: diagrammatic depiction of insole shape and internal volume of the last.


Multi-support shoes

Raised upper – plus a sturdy arthrodesis cap

The multi-support shoe has an even more raised upper than our support shoe, plus a sturdy arthrodesis cap. This area can be further reinforced through the use of mouldable thermo-plastic support strips on the inside and outside. These stiffening strips are available in two strengths: strong (white) and medium (skin-coloured). The shoe comes with four white strips. These provide even more lateral and rotational stability in the ankle joint area without hampering the natural rolling of the foot. The integrated rocker bar makes walking easier.

As with our other therapeutic shoes, individual shoe adjustments can be made with these models too. The applications for the multi-support shoe include the ongoing treatment of ankle joint capsular ligament injuries which originally received conservative or surgical treatment, as well as capsular ligament insufficiency and Achilles tendon injuries. It is also suitable for immobilising the tarsal and mid-foot regions and for the ongoing treatment of fractures and arthrodesis.





Where each foot requires different care,
shoes of differing sizes can be combined.

Anti-varus shoes

Reliable foot care right from the start

Childhood pigeon toe is a deformation which is frequently congenital and responds well to conservative treatment. Pigeon toe treatment should begin as early as possible in order to prevent long-term harm, so suitable countermeasures are taken from infancy onwards.

Depending on the severity of the misalignment, corrective treatment may be taken in the form of physiotherapy or the use of plaster casts or night splints.

If the misalignment persists until the child starts walking, corrective anti-varus shoes combined with insoles may be used.

The indication for anti-varus shoes is the treatment of walking-age children with forefoot adduction such as pigeon toe, serpentine foot or previously treated club foot.

Expectations for therapy, comfort and normal shoe appearance are met.

Successfully
guiding growth





The effectiveness of our anti-varus shoes is founded on our proven three-point correction system:

The first two correction points, on the inside of the heel and at the level of the Lisfranc joint, support the rear foot, thus protecting against torsion in the ankle joint during the necessary forefoot correction process.

The third correction point is on the inside of the big toe's metatarsophalangeal joint. This promotes correction of the forefoot and prevents the foot from reverting to the deformed position.

To bring about these corrective measures, the anti-varus shoe features a heel cap which extends along the inner side of the shoe. The corrective properties of our anti-varus shoes support successful conservative treatment of pigeon toe. Correction comes into effect where there is outward axial deviation of at least 8°.



Where each foot requires different care,
shoes of differing widths and sizes can be combined.

Insole care

Insoles are often vital when using therapeutic shoes.

Our NovaPED insole blanks are available in various designs and material combinations. Here we present a small selection from our insoles range for use with a variety of different indications.

To correct pigeon toe or previously treated club foot, a formed thermoplastic shell-shaped insole with prolonged medial part can be used.

The S90 shell-shaped insoles, with its tried-and-tested foil technique, has a supportive, strain-relieving and corrective effect. For instance, the insoles may be used for activating muscle tension in the case of slack pes valgus or to reduce the muscle tension treating a pes cavus or to improve the gait pattern in infantile cerebral palsy patients.

NovaPED also offers suitable assistance in influencing muscular tension. For example, the insoles are used to stimulate muscular tension in limp, skewed and sunken feet, to reduce muscular tension in cases of high arches or to improve the gait in patients suffering from infantile cerebral palsy.

Item no.	Design
440930-516	Shell-shaped insole, big toe loop, Supralen, Multiform cover, red
952480-712	Shell-shaped insole, Supinara, short sole, foil technique, Multiform cover, pink
470640-000	Shell-shaped insole, Peddy, short sole, Multiform cover, yellow
98 0649-000	Shell-shaped insole, long sole, foil technique, Multiform cover, orange
98 0260-000	Shell-shaped insole, long sole, heel plateau, foil technique, Multiform, orange
982260-207	Shell-shaped insole, posterisation, long sole, S90 foil technique, Multicoloured soft cushioned cover
952130-000	Shell-shaped insole, long sole, foil technique, microfibre soft cushioned cover, Multiform, blue
770340-000	Milled proprioceptive insole, long sole, made of ARUcut 25, with support sole





Item no.: 44 0930-516



Item no.: 98 0649-000



Item no.: 95 2480-712



Item no.: 98 0260-000



Item no.: 47 0640-000



Item no.: 98 2260-207



Item no.: 95 2130-000



Item no.: 77 0340-000

Quality management

Sustainable,
eco-friendly production



For many years, our company has been committed to sustainable, environmentally friendly production. This applies to sourcing the leather as much as its further processing: for a safer, cleaner environment.

TÜV-certified

As a conscientious family business, it is a matter of course that we manufacture according to exemplary quality standards. The process includes independent, external quality audits to which we voluntarily subject ourselves every year. Since 1997 the quality management certificate issued by the 'TÜV Rheinland' has attested to the efficiency of our enterprise. The certificate confirms that our working processes are systematised, thereby guaranteeing the fabrication of high-quality products.

Materials

Leather is part of nature

We work exclusively with manufacturers who meet our high expectations regarding high-quality workmanship and sustainable production. That is the only way to ensure that we fully satisfy our customers' health care needs.

The leather we use for our therapeutic shoes is therefore subject to statutory testing requirements. Permanent controls from the tanning stage to the finished shoe ensure consistently excellent quality. Because of this, minor scratches or scar tissue from insect bites, which you may find on the upper of the leather, are clear indications of its natural source rather than defects in workmanship.

You can have maximum confidence in the quality of the leather, we use thanks to the leather symbol certification mark.



Leather care

Children's shoes are often subject to extreme stresses and, in this process, the leather uppers must meet certain challenges. As a result, regular care is advisable.

A shoe polish can be used for nappa leather, whereas nubuck leather has a velvety surface which can become grubby over time. Rubbing with an abrasive pad or eraser is a good way of improving its appearance, while special nubuck pens will protect the leather upper.


Before putting them on for the first time, the upper material should be treated with a wetness protection spray to keep the dampness out. As a general rule, always dry wet shoes slowly at room temperature rather than exposing them to an additional heat source.



Materials

We've got something
against wet feet

Nothing is worse than cold, clammy feet. To guard against this, SympaTex has developed a special membrane which absorbs sweat and conveys it away from the shoes, so that your feet remain warm and dry all the time. Conversely, the SympaTex membrane prevents wetness getting in from the outside, thus providing effective protection whilst also being breathable.





Sole types

Orthotic, support, anti-varus and neutral shoes

- The flexible outsole promotes the foot's flexing action
 - In our anti-varus and neutral shoes, the outsole is less flexible due to the integrated three-point correction
 - Shell bottom with deep space to accommodate foot (approx. 15 mm at size 28)
 - Wide contact area
 - Foamed PU bottom (full material)
 - Easy to modify
-
- Applied leather frame
 - Wide contact area
 - Foamed EVA bottom in ray pattern (the medial extended heel supports the longitudinal arch)
 - Easy to modify
-
- Rocker bottom
 - Wider contact area
 - Micropore midsole
 - The triangular profile of the rubber outsole ensures a secure stance
 - Very easy to modify



Please note that the toe spring, outsole production and sole shape and width may vary, even within one shoe segment.

Support shoes

- The flexible profile sole supports a natural gait
- Shell bottom with extra depth (approx. 10 mm at size 28)
- Wider contact area
- Foamed PU bottom with air chamber system
- Modifiable



Anti-varus and neutral shoes

- Shell bottom with extra depth (approx. 10 mm at size 28)
- Wide contact area
- Foamed PU bottom with air chamber system
- Modifiable

Sole materials, heels, etc.

That's how trendy outsoles can be

Sole sheets are available in size 38 × 26 cm
and in the following colours.

Item no.	Colour
231940-000	black
231940-001	white
231940-002	honey
231940-003	pink
231940-004	blue
231940-005	orange
231940-006	mamba green
231940-007	grey
231940-008	red




Rotation heels

At every step, the (de-)rotation heel aids the natural gait and facilitates the correction of leg rotation problems.

The effective wearing time for a pair of (de-) rotation heels ranges from four to six weeks. If individual ribs break off during that period, rotation will not be diminished – even three or four ribs allow effective rotation. You can check this for yourself by manually pressing down the shoe in the heel area.

Item no.	Form	Size (approx.)	Colour
012250-010	round	68 x 68 mm	black
012250-020	round	87 x 87 mm	black

Correction generally takes around nine months, during which time at least six pairs of (de-)rotation heels will be needed. Naturally, the treatment time depends on the severity of the rotation problem.



Sole materials, heels, etc.

Shoelaces in a variety of colours and lengths

Subtle shades or bold combinations: the choice is yours. We offer shoelaces in many different colours and lengths, allowing you to individually design your shoes. Available in lengths 75, 90, 120 and 150 cm.

In our 'my design' configurator you can view the various shoelaces first hand in different shoes and combine them yourself.
www.mydesign-schein.de



Item no.	Colour
239000	black
239001	white
239002	dark blue
239003	orange
239004	mamba green
239005	royal
239006	pink
239007	red
239008	anthracite
239009	grey
239010	grey glitter
239011	white glitter
239012	blue glitter
239013	bordeaux glitter

The Schein app

The safe way to the right shoe

The Schein app - our quick finder service

From sizing to the selection of all available shoe models, our app supports you in providing advice every single day.

You have the opportunity to choose between specifying the last or ascertaining warehouse availability. This makes it even easier to check the availability of our my generation therapeutic shoes. Simply enter the model number you are looking for, select the desired width or size and our system will immediately show you the warehouse availability.


Your benefits at a glance:

- No waiting time on the phone
- Order processing in just a few clicks
- Simplified workflow

Registration:



Select treatment:



Width selection:




Search results:



Stock availability:



Online ordering:



my design configurator

My individual shoe

The design of our therapeutic shoes is strongly influenced by everyday shoes, thus allowing us to incorporate the fashionable components which are so important, particularly among younger people. Thanks to the diverse range of materials, leather types and materials we use, we have created an extensive collection which combines fashion with therapeutic benefits to provide feet with ideal care.

However, we do not always satisfy everybody's personal taste. For instance, children, teenagers and adults like a particular model, but would rather have it in a combination of colour or material different to the one in our catalogue.

That is why we have developed the 'my design' shoe configurator, which allows you to create your own unique pair of shoes entirely to suit your own personal tastes.

From the sole to the shoelaces, there are many different design variations, plus a wide range of colours. Just print out your chosen design and take it to your specialised dealer to order the model you have configured. We reserve the right to make changes to colours and designs.

You can recognise the selected models, which you have configured according to your wishes, by looking for the following symbol:



The screenshot shows the 'my design' shoe configurator interface. At the top, there's a navigation bar with 'myGeneration' on the left, the 'my design' logo in the center, and 'Ultimate shoes', 'Support shoes', and 'Anti-varus shoes' on the right. Below the logo is a small 'deutsch' button. The main area has tabs for '1. Treatment', '2. Model', '3. Size / Width', and '4. Shoe Design'. The 'Shoe Design' tab is active, showing a 3D rendering of a blue and white high-top sneaker with three straps. To the right of the shoe is a legend titled 'Surfaces' listing various parts: Hook-and-loop fasteners (light grey), Surface 1 (dark grey), Surface 2 (medium grey), Tongue (dark grey), Border (dark grey), Cap (light grey), Emblem (dark grey), Eyelets (light grey), Inner lining (dark red), Seams (dark red), and Sole (light grey). At the bottom left, there's a 'CUSTOMER or RETAILER' section with a note about saving configurations. At the bottom right are buttons for 'Save', 'Continue', and 'Order'.

my design models



Orthotic shoe: catalogue pages 48–49
 Support shoe: catalogue pages 97–98
 Anti-varus shoe: catalogue pages 123–124



	Item no.	Width	Size	Fastener	Reflector
Orthotic shoe	379300	TN8	18–42	Hook-and-loop fasteners with roller loop	yes
Support shoe	357000	TN4	23–42	Hook-and-loop fasteners with roller loop	yes
Anti-varus shoe	341000 341005	AV NS	18–42	Hook-and-loop fasteners with roller loop	yes

my design models



Orthotic shoe: catalogue pages 44–45
Anti-varus shoe: catalogue pages 119

Configuration examples
We reserve the right to make changes to colours and designs



Orthotic shoe: catalogue page 60

Configuration examples
We reserve the right to make changes to colours and designs

	Item no.	Width	Size	Fastener	Reflector
Orthotic shoe	361100 371100	TN5 TN8	21–45	Hook-and-loop fasteners with no roller loop	yes
Anti-varus shoe	347000 347005	AV NS	21–35	Hook-and-loop fasteners with no roller loop	yes
Orthotic shoe	366500 376500	TN5 TN8	20–47	Lacing	yes



Configuration examples
We reserve the right to make changes to colours and designs



	Item no.	Width	Size	Fastener	Reflector
Orthotic shoe	365000	TN5	20–50	Hook-and-loop fasteners with roller loop	yes
	375000	TN8			
	375005	TN9			
	375055	TN12			
Orthotic shoe	368300	TN5	36–47	Lacing	
	378300	TN8			yes

my design models



Configuration examples
We reserve the right to make changes to colours and designs



Configuration examples
We reserve the right to make changes to colours and designs

	Item no.	Width	Size	Fastener	Reflector
Orthotic shoe	368390 378390	TN5 TN8	24–45	Lacing / Zip on the inside and outside	no
Orthotic shoe	368470 378470	TN5 TN8	21–47	Lacing / Hook-and-loop fasteners	no



Orthotic shoe: catalogue pages 78–79
Support shoe: catalogue pages 102 / 103
Anti-varus shoe: catalogue page 133



Configuration examples
We reserve the right to make changes to colours and designs



Orthotic shoe: catalogue page 68
Support shoe: catalogue pages 100–101
Anti-varus shoe: catalogue pages 126–128



Configuration examples
We reserve the right to make changes to colours and designs

	Item no.	Width	Size	Fastener	Reflector
Orthotic shoe	362200 372200	TN5 TN8	20–47	Hook-and-loop fasteners, top one with roller loop	yes
Support shoe	356000	TN2	23–47	Hook-and-loop fasteners, top one with roller loop	yes
Anti-varus shoe	342000 342005	AV NS	20–42	Hook-and-loop fasteners, top one with roller loop	yes
Orthotic shoe	368000 378000	TN5 TN8	20–47	Lacing / Zip on the inside	yes
Support shoe	356100	TN2	23–47	Lacing / Zip on the inside	yes
Anti-varus shoe	343000 343005	AV NS	20–42	Lacing	yes

Orthotic shoes

Our orthotic shoes – a combination of proven medical functionality and a fresh appearance

The right footwear is an important addition to an individually produced orthosis. Shoes worn by people wearing an orthosis not only have to satisfy certain key requirements, but must also make walking and running easier.

For these reasons, we remain in close contact with our customers with the aim of making continuous improvements and refinements:

- Orthotic shoe width TN12
- Orthotic shoes for joint orthoses
- Flexible outsole
- New rear access
- Straight posture

Due to their material quality and shape, conventionally manufactured shoes are not up to meeting the requirements of orthoses. That is why we recommend our therapeutic shoes.



Orthotic shoe width TN12

To facilitate treatment involving bulky orthoses, we have developed – together with our craftsmen – the width TN12 for our orthotic shoes.



We supply every shoe in the width TN12 with two insoles in 40 and 50 Shore, meaning that the volume can be changed by removing one or both of the foot orthotics.

This ensures that the external design of the shoe is not affected. If a therapeutic shoe is only needed for one foot, a different width and shoe size can be provided for the normal foot, in the same way as all the other TN widths in the orthotic shoe category.



Orthotic shoes for joint orthoses

Wearers of joint orthoses need a shoe with a specially shaped heel. Although the shoe will have a more voluminous rear-foot and heel area, it must nevertheless fit well. That is the only way to ensure it holds the foot properly in place and guarantees correct foot movement when walking.

Our orthotic shoes, which are marked with the symbol shown on the right, have been specially designed to meet these requirements.



Flexible outsole

Ever thinner and more flexible materials are used in the orthopaedic handcraft and its treatments. In response, we have introduced an even more flexible version of the sole for our orthotic and support shoes.

All shoes designated in the tables as 'flexible' are extremely elastic. None of this has any negative impact on the therapeutic benefits of the orthosis.

Orthotic shoes

Sandals

Page 40–49



Low shoes

Page 50–67



Boots

Page 68–85



School shoes

Page 87



Limited quantities

Page 88–91



Orthotic shoes



Orthotic shoes



Opens
very wide

Pippa

367020 / width TN5
377020 / width TN8



Samara

367060 / width TN5
377060 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Pippa	367020 377020	TN5 TN8	21 - 32	Nubuck	turquoise	Hook-and-loop fasteners with roller loop	textile	flexible	yes
Samara	367060 377060	TN5 TN8	21 - 32	Nubuck	shiny silver	Hook-and-loop fasteners with roller loop	textile	flexible	yes



Orthotic shoes

Orthotic shoes



Chen

367000 / width TN5
377000 / width TN8
377005 / width TN12



Niklas

367080 / width TN5
377080 / width TN8

Opens
very wide



Paula

367030 / width TN5
377030 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Chen	367000	TN5						
	377000	TN8	20 - 45	Nubuck	dark blue / black	Hook-and-loop fasteners with roller loop	textile	yes
	377005	TN12						
Paula	367030	TN5	20 - 38	Nubuck / Mesh fabric	grey / pink	Hook-and-loop fasteners with roller loop	textile	yes
	377030	TN8						
Niklas	367080	TN5	20 - 45	Nubuck / Mesh fabric	light grey / dark blue	Hook-and-loop fasteners with roller loop	textile	yes
	377080	TN8						



Camillo

367050 / width TN5
377050 / width TN8



Nell

367070 / width TN5
377070 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Camillo	367050 377050	TN5 TN8	24 - 35	Nappa / Faded denim look	blue	Hook-and-loop fasteners with roller loop	textile	flexible	yes
Nell	367070 377070	TN5 TN8	24 - 35	Nubuck / Fabric	dark grey / green	Hook-and-loop fasteners with roller loop	textile	flexible	yes

Orthotic shoes



You can use the configurator to design this model for yourself



Carbon

361900 / width TN5
371900 / width TN8



Marina

361910 / width TN5
371910 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Carbon	361900 371900	TN5 TN8	21 - 45	Nappa / Mesh fabric	carbon black	Hook-and-loop fasteners with no roller loop	calf leather lining	yes
Marina	361910 371910	TN5 TN8	21 - 32	Nubuck	turquoise	Hook-and-loop fasteners with no roller loop	calf leather lining	yes



Daisy

361140 / width TN5

Mareike

363960 / width TN5

373960 / width TN8



Kletty

363950 / width TN5

373950 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Daisy	361140	TN5	21 - 28	Nappa	black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Kletty	363950 373950	TN5 TN8	21 - 35	Nubuck / Mesh fabric	dark blue	Hook-and-loop fasteners with no roller loop	calf leather lining	yes
Mareike	363960 373960	TN5 TN8	21 - 32	Nubuck / Mesh fabric	pink	Hook-and-loop fasteners with no roller loop	calf leather lining	yes

Orthotic shoes



Grey

371880 / width TN8



Anouk

371860 / width TN8



Sofia

371830 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Sofia	371830	TN8	21 - 32	Patent leather embossed	red	Hook-and-loop fasteners with roller loop	calf leather lining	no
Anouk	371860	TN8	21 - 32	Nubuck / Nappa	shiny silver / white	Hook-and-loop fasteners with roller loop	calf leather lining	no
Grey	371880	TN8	21 - 35	Nubuck / Nappa	grey	Hook-and-loop fasteners with roller loop	calf leather lining	no



Orthotic shoes

Orthotic shoes



Hero

361080 / width TN5
371080 / width TN8



Domino

371050 / width TN8



Emilia

363040 / width TN5
373040 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Hero	361080	TN5	18 - 45	Nubuck / Mesh fabric	dark blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes
	371080	TN8	18 - 45					
Emilia	363040	TN5	21 - 35	Nubuck	pink	Hook-and-loop fasteners with roller loop	calf leather lining	yes
	373040	TN8						
Domino	371050	TN8	27 - 42	Nappa / Nubuck	grey	Hook-and-loop fasteners with roller loop	calf leather lining	yes



Jack

361020 / width TN5
371020 / width TN8

Mila

371070 / width TN8



You can use the configurator to design this model for yourself



my design

www.mydesign-schein.de

Tuffy

371090 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Jack	361020 371020	TN5 TN8	21 - 38	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Mila	371070	TN8	21 - 35	Patent leather embossed	red	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Tuffy	371090	TN8	24 - 42	Nubuck / Nappa	grey / blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes

Orthotic shoes



Klara

362570 / width TN5

372570 / width TN8



Ocean

362550 / width TN5

372550 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Ocean	362550 372550	TN5 TN8	18 - 35	Faded denim look / Nubuck	blue / brown	Hook-and-loop fastener with roller loop	calf leather lining	flexible	no
Klara	362570 372570	TN5 TN8	18 - 32	Faded denim look / Nubuck	pink	Hook-and-loop fastener with roller loop	calf leather lining	flexible	no

A close-up, profile photograph of a young boy with light brown hair, lying on his stomach on a light-colored, textured carpet. He is wearing a white t-shirt and has his head propped up by his hands. He is looking upwards and slightly to the left. The background is blurred, showing what appears to be a window or a doorway.

Orthotic shoes

Orthotic shoes



Amigo

364220 / width TN5
374220 / width TN8

Kate

364290 / width TN5
374290 / width TN8



Navaro

364200 / width TN5
374200 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Navaro	364200 374200	TN5 TN8	30 - 47	Nappa	black	Hook-and-loop fasteners with no roller loop	textile	flexible	no
Amigo	364220 374220	TN5 TN8	24 - 35	Faded denim look	blue	Hook-and-loop fasteners with no roller loop	textile	flexible	no
Kate	364290 374290	TN5 TN8	24 - 32	Nubuck with floral embossing	shiny silver	Hook-and-loop fasteners with no roller loop	textile	flexible	no



Rio

364210 / width TN5
374210 / width TN8

Nora

364270 / width TN5
374270 / width TN8



Bailey

364260 / width TN5
374260 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Rio	364210 374210	TN5 TN8	24 - 39	Nubuck / Nappa	blue / dark blue	Hook-and-loop fasteners with roller loop	textile	flexible	no
Bailey	364260 374260	TN5 TN8	30 - 42	Nappa	white	Hook-and-loop fasteners with roller loop	textile	flexible	no
Nora	364270 374270	TN5 TN8	24 - 32	Nubuck	turquoise	Hook-and-loop fasteners with roller loop	textile	flexible	no

Orthotic shoes



Fox

366780 / width TN5
376780 / width TN8



Ella

366770 / width TN5
376770 / width TN8



Finja

366760 / width TN5
376760 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Finja	366760 376760	TN5 TN8	21 - 32	Nubuck	shiny silver	Hook-and-loop fasteners with roller loop	textile	yes
Ella	366770 376770	TN5 TN8	21 - 35	Nubuck	lilac	Hook-and-loop fasteners with roller loop	textile	yes
Fox	366780 376780	TN5 TN8	21 - 47	Nubuck	grey	Hook-and-loop fasteners with roller loop	textile	yes



Orthotic shoes

Orthotic shoes



Tim

366750 / width TN5
376750 / width TN8



Linda

366740 / width TN5
376740 / width TN8



Paco

366720 / width TN5
376720 / width TN8
376725 / width TN9
378725 / width TN12



Van

368700 / width TN5
378700 / width TN8

You can use the configurator to design this model for yourself



	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Paco	366720	TN5						
	376720	TN8						
	376725	TN9	20 - 47	Nubuck / Mesh fabric	grey / dark blue	Hook-and-loop fasteners with roller loop	textile	yes
	378725	TN12						
Linda	366740	TN5						
	376740	TN8	20 - 32	Nubuck / Mesh fabric	grey / pink	Hook-and-loop fasteners with roller loop	textile	yes
Tim	366750	TN5						
	376750	TN8	20 - 35	Nubuck / Mesh fabric	grey / mint	Hook-and-loop fasteners with roller loop	textile	yes
Van	368700	TN5						
	378700	TN8	36 - 50	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	textile	yes



Blue

363320 / width TN5
373320 / width TN8

Orthotic shoes

Maya

363310 / width TN5
373310 / width TN8



Buddy

363750 / width TN5
373750 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Maya	363310 373310	TN5 TN8	21 - 38	Nubuck / Mesh fabric	pink	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Blue	363320 373320	TN5 TN8	21 - 42	Nubuck / Mesh fabric	dark blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Buddy	363750 373750	TN5 TN8	26 - 42	Nubuck / Mesh fabric	dark blue / grey	Hook-and-loop fasteners with no roller loop	textile	yes

Orthotic shoes



Arne

363760 / width TN5
373760 / width TN8



Fynn

362810 / width TN5
370810 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Fynn	362810 370810	TN5 TN8	30 - 47	Nappa	black	Hook-and-loop fasteners with no roller loop	textile	yes
Arne	363760 373760	TN5 TN8	26 - 42	Nappa / Mesh fabric	black / grey	Hook-and-loop fasteners with no roller loop	textile	yes



Blacky

363030 / width TN5
373030 / width TN8



Gym

363010 / width TN5
373010 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Gym	363010 373010	TN5 TN8	26 - 47	Nappa / Mesh fabric	white	Lacing	calf leather lining	yes
Blacky	363030 373030	TN5 TN8	30 - 47	Nappa / Nubuck / Mesh fabric	black	Lacing	calf leather lining	yes

Orthotic shoes



Mathilda

366660 / width TN5

376660 / width TN8



Nino

366620 / width TN5

376620 / width TN8



Indigo

366640 / width TN5

376640 / width TN8

You can use the
configurator to design
this model for yourself



	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Nino	366620 376620	TN5 TN8	20 - 47	Nubuck / Mesh fabric	grey / dark blue	Lacing	textile	yes
Indigo	366640 376640	TN5 TN8	21 - 47	Nubuck	dark blue	Lacing	textile	yes
Mathilda	366660 376660	TN5 TN8	24 - 35	Nubuck	shiny silver	Lacing	textile	yes

A portrait of a young man with dark hair, smiling broadly with his arms crossed. He is wearing a white t-shirt with a black bowtie graphic and a blue denim jacket over it. A pair of black headphones hangs around his neck. The background is a solid yellow color.

Orthotic shoes

Orthotic shoes



Neo

369660 / width TN5
379660 / width TN8



Sam

369620 / width TN5
379620 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Sam	369620 379620	TN5 TN8	36 - 45	Nappa / Nubuck	dark blue / grey	Lacing	textile	flexible	no
Neo	369660 379660	TN5 TN8	30 - 42	Nappa	white	Lacing	textile	flexible	no



Anton

369680 / width TN5

379680 / width TN8



Theo

369600 / width TN5

379600 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Theo	369600 379600	TN5 TN8	30 - 45	Nubuck	black	Lacing	textile	flexible	no
Anton	369680 379680	TN5 TN8	30 - 45	Nubuck	grey / blue	Lacing	textile	flexible	no

Orthotic shoes



Noel

369780 / width TN5
379780 / width TN8

Lotta

369710 / width TN5
379710 / width TN8



Adrian

369700 / width TN5
379700 / width TN8



Phil

369720 / width TN5
379720 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Adrian	369700 379700	TN5 TN8	24 - 39	Nubuck	black	Lacing / Zip on the inside and outside	textile	flexible	no
Lotta	369710 379710	TN5 TN8	21 - 32	Nubuck	lilac	Lacing / Zip on the inside and outside	textile	flexible	no
Phil	369720 379720	TN5 TN8	24 - 32	Nubuck	blue	Lacing / Zip on the inside and outside	textile	flexible	no
Noel	369780 379780	TN5 TN8	24 - 39	Nubuck / Nappa	grey	Lacing / Zip on the inside and outside	textile	flexible	no




Orthotic shoes

Orthotic shoes



Nanaimo

368220 / width TN5
378220 / width TN8



Trivor

368200 / width TN5
378200 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Trivor	368200 378200	TN5 TN8	36 - 45	Nappa	black	Lacing / Zip on the inside and outside	textile	flexible	no
Nanaimo	368220 378220	TN5 TN8	24 - 42	Faded denim look	blue	Lacing / Zip on the inside and outside	textile	flexible	no



You can use the configurator to design this model for yourself



Fabienne

368290 / width TN5
378290 / width TN8

Alina

368280 / width TN5
378280 / width TN8

Orthotic shoes

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Alina	368280 378280	TN5 TN8	24 - 39	Faded denim look	grey	Lacing / Zip on the inside and outside	textile	flexible	no
Fabienne	368290 378290	TN5 TN8	24 - 32	Faded denim look	lilac	Lacing / Zip on the inside and outside	textile	flexible	no

Orthotic shoes

Applies to all models
on this double page



Lina

366360 / width TN5
376360 / width TN8



Yannie

364300 / width TN5
374300 / width TN8



Matti

366390 / width TN5
376390 / width TN8




You can use the
configurator to design
this model for yourself



	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Yannie	364300 374300	TN5 TN8	20 - 47	Faded denim look	black	Lacing / Zip on the inside	textile	—	yes
Lina	366360 376360	TN5 TN8	24 - 42	Nubuck with floral embossing	shiny silver	Lacing / Zip on the inside	textile	flexible	yes
Matti	366390 376390	TN5 TN8	30 - 47	Nappa	black	Lacing / Zip on the inside	textile	flexible	yes

- Easy access through the movable, fold-out back cap
- Adjustment to suit foot possible

Applies to all models on this page



You can use the configurator to design this model for yourself



Milow

369200 / width TN5
379200 / width TN8



Levi


369260 / width TN5
379260 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Milow	369200 379200	TN5 TN8	21 - 47	Faded denim look	black	Lacing / Hook-and-loop fastener	textile	flexible	no
Levi	369260 379260	TN5 TN8	21 - 47	Faded denim look	white	Lacing / Hook-and-loop fastener	textile	flexible	no

Orthotic shoes



Applies to all models on this double page



Himba

369240 / width TN5
379240 / width TN8



Vio

369220 / width TN5
379220 / width TN8

Elise

369230 / width TN5
379230 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Vio	369220	TN5	21 - 47	Faded denim look	blue	Lacing / Hook-and-loop fastener	textile	flexible	no
	379220	TN8							
Elise	369230	TN5	21 - 35	Faded denim look	pink	Lacing / Hook-and-loop fastener	textile	flexible	no
	379230	TN8							
Himba	369240	TN5	21 - 35	Faded denim look	yellow	Lacing / Hook-and-loop fastener	textile	flexible	no
	379240	TN8							

- Easy access through the movable, fold-out back cap
- Adjustment to suit foot possible



Nala

369250 / width TN5
379250 / width TN8



You can use the configurator to design this model for yourself



Malia

369210 / width TN5
379210 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Malia	369210 379210	TN5 TN8	21 - 39	Nubuck	lilac	Lacing / Hook-and-loop fastener	textile	flexible	no
Nala	369250 379250	TN5 TN8	21 - 35	Nubuck	turquoise	Lacing / Hook-and-loop fastener	textile	flexible	no

Orthotic shoes



Lynn


362030 / width TN5



Louis

360700 / width TN5

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Louis	360700	TN5	23 - 42	Nappa	black	Lacing	calf leather lining	no
Lynn	362030	TN5	21 - 35	Nappa / Nubuck	black	Lacing	calf leather lining	no

A close-up portrait of a young boy with short brown hair, smiling broadly at the camera. He is wearing a dark blue school blazer over a white collared shirt and a red, white, and black striped tie. The background is blurred, showing what appears to be a school hallway or entrance. The lighting is warm and focused on his face.

Orthotic shoes

Orthotic shoes



Ole

364920 / width TN5
374920 / width TN8



Luana

364910 / width TN5
374910 / width TN8



Lumi

364940 / width TN5
374940 / width TN8


	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Luana	364910 374910	TN5 TN8	21 - 35	Nubuck	lilac	Hook-and-loop fasteners, lower one with no roller loop	textile	flexible	no
Ole	364920 374920	TN5 TN8	21 - 32	Nubuck	royal blue	Hook-and-loop fasteners, lower one with no roller loop	textile	flexible	no
Lumi	364940 374940	TN5 TN8	21 - 35	Nubuck	ocker	Hook-and-loop fasteners, lower one with no roller loop	textile	flexible	no



Ava

364970 / width TN5
374970 / width TN8

Applies to all models
on this double page



Ronja

364950 / width TN5
374950 / width TN8




River

364980 / width TN5
374980 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Ronja	364950 374950	TN5 TN8	21 - 29	Nappa	red	Hook-and-loop fasteners, lower one with no roller loop	textile	flexible	no
Ava	364970 374970	TN5 TN8	21 - 35	Nubuck	turquoise	Hook-and-loop fasteners, lower one with no roller loop	textile	flexible	no
River	364980 374980	TN5 TN8	30 - 47	Nubuck / Nappa	grey / dark grey	Hook-and-loop fasteners, lower one with no roller loop	textile	flexible	no

Orthotic shoes

Applies to all models
on this page



Straight posture



Glitzi

362660 / width TN5
372660 / width TN8



Unique

362680 / width TN5
372680 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Glitzi	362660 372660	TN5 TN8	24 - 32	Nubuck	shiny silver	Hook-and-loop fasteners, top one with roller loop	textile	flexible	no
Unique	362680 372680	TN5 TN8	27 - 42	Nubuck / Nappa	black / grey	Hook-and-loop fasteners, top one with roller loop	textile	flexible	no

A vibrant, candid photograph of five young women of diverse ethnicities, all smiling broadly at the camera. They are posed in a close-knit, circular arrangement, with some individuals leaning their heads against each other. The woman in the foreground is wearing a white and blue horizontally striped t-shirt. The background is a bright, overexposed outdoor scene with a person's legs and jeans visible.

Orthotic shoes

Orthotic shoes



You can use the configurator to design this model for yourself



my design
www.mydesign-schein.de

Applies to all models on this page



Curly

364110 / width TN5
374110 / width TN8



Julian

364170 / width TN5
374170 / width TN8
374175 / width TN12

Maris

364140 / width TN5
374140 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Curly	364110 374110	TN5 TN8	20 - 42	Faded denim look	blue	Hook-and-loop fasteners, top one with roller loop	textile	yes
Maris	364140 374140	TN5 TN8	23 - 45	Nappa	brown	Hook-and-loop fasteners, top one with roller loop	textile	yes
Julian	364170 374170 374175	TN5 TN8 TN12	20 - 42	Nappa	dark blue	Hook-and-loop fasteners, top one with roller loop	textile	yes



Toby

360770 / width TN5
370770 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Toby	360770 370770	TN5 TN8	23 - 38	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Jacko	364180 374180	TN5 TN8	23 - 38	Nappa	black	Hook-and-loop fasteners, top one with roller loop	calf leather lining	yes

Orthotic shoes



Dani

363460 / width TN5
373460 / width TN8



Ida

363430 / width TN5
373430 / width TN8



Samu

363410 / width TN5
373410 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Samu	363410 373410	TN5 TN8	21 - 45	Faded denim look / Nappa	grey	Hook-and-loop fasteners with roller loop	textile	yes
Ida	363430 373430	TN5 TN8	21 - 32	Faded denim look	red	Hook-and-loop fasteners with roller loop	textile	yes
Dani	363460 373460	TN5 TN8	24 - 42	Nubuck	blue	Hook-and-loop fasteners with roller loop	textile	yes



Vegas

363420 / width TN5
373420 / width TN8



Stella

363450 / width TN5
373450 / width TN8
373455 / width TN9
378455 / width TN12



Phoenix

363480 / width TN5
373480 / width TN8
373485 / width TN9
378485 / width TN12

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Vegas	363420 373420	TN5 TN8	24 - 45	Nubuck / Mesh fabric	dark blue / grey	Hook-and-loop fasteners with roller loop	textile	yes
Stella	363450 373450 373455 378455	TN5 TN8 TN9 TN12	21 - 38	Nubuck	purple / grey	Hook-and-loop fasteners with roller loop	textile	yes
Phoenix	363480 373480 373485 378485	TN5 TN8 TN9 TN12	18 - 50	Nubuck / Mesh fabric	black / grey	Hook-and-loop fasteners with roller loop	textile	yes

Orthotic shoes



Micki

363530 / width TN5



Robin

363510 / width TN5
373510 / width TN8



Matterhorn

363620 / width TN5

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Robin	363510 373510	TN5 TN8	36 - 47	Nappa / Nubuck / Mesh fabric	black / grey / petrol	Lacing / Hook-and-loop fastener with roller loop	calf leather lining	yes
Micki	363530	TN5	36 - 47	Nappa / Mesh fabric	brown	Lacing / Hook-and-loop fastener with roller loop	calf leather lining	yes
Matterhorn	363620	TN5	30 - 47	Nappa / Mesh fabric	black	Hook-and-loop fasteners, central one with no roller loop	calf leather lining	yes



Mats

362700 / width TN5
372700 / width TN8



Charlie

362720 / width TN5
372720 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Mats	362700 372700	TN5 TN8	36 - 45	Nappa	carbon black	Lacing / Zip on the inside	teddy fur	flexible	yes
Charlie	362720 372720	TN5 TN8	27 - 42	Faded denim look	blue	Lacing / Zip on the inside	teddy fur	flexible	yes

Orthotic shoes



Nevis

362600 / width TN5
372600 / width TN8



Colima

362630 / width TN5
372630 / width TN8

Mountain

362610 / width TN5
372610 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Sole	Reflector
Nevis	362600 372600	TN5 TN8	20 - 47	Faded-look waxed denim	black	Hook-and-loop fasteners, top one with roller loop	teddy fur	—	yes
Mountain	362610 372610	TN5 TN8	36 - 45	Nubuck	brown	Hook-and-loop fasteners with roller loop	teddy fur	flexible	yes
Colima	362630 372630	TN5 TN8	24 - 38	Nubuck	purple	Hook-and-loop fasteners with roller loop	teddy fur	flexible	yes



Berni

372510 / width TN8



Felix

362500 / width TN5
372500 / width TN8

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Felix	362500 372500	TN5 TN8	26 - 42	Nubuck / Mesh fabric	dark blue / grey	Hook-and-loop fasteners with roller loop	teddy fur	no
Berni	372510	TN8	26 - 38	Nubuck / Mesh fabric	red / dark blue	Hook-and-loop fasteners with roller loop	teddy fur	no

The perfect twist!

Orthotic shoes with Boa® fastening system

One twist for
good support



You can also request our separate brochure:

my generation
by scholl

Orthotic shoes with Boa® fastening system

Our therapeutic shoes are specially designed for people who have been wearing orthotics for a long time. Our special orthopedic knowledge, modern leather and high-quality materials make for long-lasting comfort.

The Boa® fastening system follows the philosophy of our orthotic shoes thanks to its robust workmanship.

Boa® fastening system

We offer fastening with a twist!

We have fitted all the orthotic shoe models in this brochure with a Boa® fastening system. Originally developed for cyclists, it offers quick, easy, secure fits for users of orthotics with limited motor function and reduced padding. Taking off the Boa® fastening system is just as easy as putting it on. This is how the Boa® fastening and unfastening is done with a single movement of the adjustment wheel. Turning the wheel allows the shoe to be adjusted easily and individually for the patient and goes to the millimetre.

The Boa® fastening system is a single hold and the highest level of comfort when wearing the shoe.

The secure hold

We call the Boa® fastening system a 'one-twist' hold. We have fitted all the orthotic shoe models in this brochure with a Boa® fastening system. Originally developed for cyclists, it offers quick, easy, secure fits for users of orthotics with limited motor function and reduced padding. Taking off the Boa® fastening system is just as easy as putting it on. This is how the Boa® fastening and unfastening is done with a single movement of the adjustment wheel. Turning the wheel allows the shoe to be adjusted easily and individually for the patient and goes to the millimetre.

Unlike conventional shoelaces, which can absorb water or grit, the low-friction synthetic coverings attached to the Boa® fastening system remain dry and clean.

Even under heavy strain, the cloth and webbing are a loop-lasting concept.

Once it has been fastened, the facing stays in place – the orthosis and/or foot are securely fixed.

How to use the Boa® fastening system:

- To fasten: just pull out the adjustment wheel and the fastening is secured.
- To unfasten: just pull out the adjustment wheel and the fastening is released.
- To adjust: push down the adjustment wheel and then set the fastening by turning the wheel clockwise.

Orthotic shoes

School shoes



Daisy

361140 / width TN5 / size 21 - 28
Page 45



Jack

361020 / width TN5 / size 21 - 38
371020 / width TN8 / size 21 - 38
Page 49



Navaro

364200 / width TN5 / size 30 - 47
374200 / width TN8 / size 30 - 47
Page 52



Van

368700 / width TN5 / size 36 - 50
378700 / width TN8 / size 36 - 50
Page 56



Fynn

362810 / width TN5 / size 30 - 47
370810 / width TN8 / size 30 - 47
Page 58



Blacky

363030 / width TN5 / size 30 - 47
373030 / width TN8 / size 30 - 47
Page 59



Trivor

368200 / width TN5 / size 36 - 45
378200 / width TN8 / size 36 - 45
Page 66



Lynn

362030 / width TN5 / size 21 - 35
Page 72



Louis

360700 / width TN5 / size 23 - 42
Page 72



Jacko

364180 / width TN5 / size 23 - 38
374180 / width TN8 / size 23 - 38
Page 79



Toby

360770 / width TN5 / size 23 - 38
370770 / width TN8 / size 23 - 38
Page 79



Matterhorn

363620 / width TN5 / size 30 - 47
Page 82

Orthotic shoes Limited quantities

Limited stocks of these models are available.

Please call us on +49 (0)2191 910 201 to ensure their availability.

Our advice: look for a suitable alternative model from our regular range.



Akado

367040 / width TN5 / size 36 - 41
377040 / width TN8 / size 40 - 44



Hannah

367090 / width TN5 / size 20 - 35
377090 / width TN8 / size 20 - 35



Jakob

367010 / width TN5 / size 23 - 34
377010 / width TN8 / size 23 - 37



Christin

371730 / width TN8 / size 20 - 30



Nicky

371700 / width TN8 / size 20 - 25



Amelie

361940 / width TN5 / size 21 - 35
371940 / width TN8 / size 21 - 35



Tolouse

361920 / width TN5 / size 21 - 35
371920 / width TN8 / size 21 - 35



Zoe

361930 / width TN5 / size 21 - 34
371930 / width TN8 / size 21 - 35



Chico

361990 / width TN5 / size 27 - 42
371990 / width TN8 / size 27 - 42



Lupi

370380 / width TN8 / size 21 - 32



Amy

370390 / width TN8 / size 27 - 31



Jessica

360910 / width TN5 / size 29 - 32
370910 / width TN8 / size 21 - 38



Tom

360900 / width TN5 / size 21 - 47
370900 / width TN8 / size 21 - 34



Pina

371060 / width TN8 / size 36 - 42



Linus

371010 / width TN8 / size 24 - 38



Lilly

371030 / width TN8 / size 30 - 35



Ben

370210 / width TN8 / size 18 - 34



Sunny

370450 / width TN8 / size 23 - 34



Sky

370520 / width TN8 / size 23 - 38



David

044165 / width TN8 / size 18 - 31



Romeo

364240 / width TN5 / size 27 - 35
374240 / width TN8 / size 26 - 32



Leonie

364230 / width TN5 / size 25 - 32
374230 / width TN8 / size 24 - 31



Strike

366710 / width TN5 / size 20 - 35
376710 / width TN8 / size 20 - 35



Edi

366730 / width TN5 / size 20 - 35



Sporty

363940 / width TN5 / size 26 - 34
373940 / width TN8 / size 27 - 33



Max

363250 / width TN5 / size 27 - 36
373250 / width TN8 / size 20 - 32



Martin

363340 / width TN5 / size 21 - 34
373340 / width TN8 / size 21 - 29



Jackson

362800 / width TN5 / size 30 - 42



Rico

362100 / width TN5 / size 36 - 42
372100 / width TN8 / size 36 - 42



Scoop

362110 / width TN5 / size 36 - 42
372110 / width TN8 / size 36 - 42



Mambo

362130 / width TN5 / size 36 - 42
372130 / width TN8 / size 36 - 42



Merle

376630 / width TN8 / size 20 - 33



Merito

368710 / width TN5 / size 36 - 47
378710 / width TN8 / size 33 - 47



Ruby

366610 / width TN5 / size 36 - 42
376610 / width TN8 / size 36 - 42



Beat

366650 / width TN5 / size 36 - 42
376650 / width TN8 / size 41 - 42



Katie

363130 / width TN5 / size 25 - 28



Olivia

363730 / width TN5 / size 26 - 28



Milan

044860 / width TN9 / size 20 - 26



Chip

363830 / width TN5 / size 30 - 41
373830 / width TN8 / size 30 - 42



Summi

363840 / width TN5 / size 30 - 38
373840 / width TN8 / size 30 - 32



Chap

363850 / width TN5 / size 30 - 42
373850 / width TN8 / size 30 - 42



Moritz

363980 / width TN5 / size 33 - 39



Jona

364420 / width TN5 / size 20 - 32
374420 / width TN8 / size 20 - 32



Cato

364440 / width TN5 / size 20 - 32
374440 / width TN8 / size 20 - 32

Orthotic shoes Limited quantities

Please call us on +49 (0)2191 910 201 to ensure their availability.



Joschi

366310 / width TN5 / size 24 - 32
376310 / width TN8 / size 24 - 32



Banu

364350 / width TN5 / size 20 - 32
374350 / width TN8 / size 20 - 30



Lara

364330 / width TN5 / size 20 - 34
374330 / width TN8 / size 22 - 32



Nicolai

366340 / width TN5 / size 30 - 42
376340 / width TN8 / size 30 - 42



Happy

366320 / width TN5 / size 24 - 38
376320 / width TN8 / size 24 - 38



Racing

366330 / width TN5 / size 24 - 32
376330 / width TN8 / size 24 - 31



Jazz

366380 / width TN5 / size 24 - 32
376380 / width TN8 / size 24 - 32



Nolly

366350 / width TN5 / size 24 - 42
376350 / width TN8 / size 24 - 42



Luisa

366370 / width TN5 / size 24 - 32
376370 / width TN8 / size 24 - 32



Filou

368110 / width TN5 / size 21 - 32
378110 / width TN8 / size 20 - 32



Jasmin

368130 / width TN5 / size 21 - 31
378130 / width TN8 / size 20 - 30



Jarvis

368100 / width TN5 / size 24 - 42
378100 / width TN8 / size 24 - 42



Momo

368120 / width TN5 / size 24 - 38
378120 / width TN8 / size 26 - 38



Weltmeister

368160 / width TN5 / size 20 - 42
378160 / width TN8 / size 20 - 42



Henry

364120 / width TN5 / size 26 - 45
374120 / width TN8 / size 26 - 31



Tina

364010 / width TN5 / size 27 - 31
374010 / width TN8 / size 27 - 29



Tommy

364000 / width TN5 / size 26 - 33
374000 / width TN8 / size 26 - 32



Elena

364930 / width TN5 / size 24 - 32
374930 / width TN8 / size 24 - 32



Mackay

364900 / width TN5 / size 40 - 47
374900 / width TN8 / size 36 - 47



Oskar

363470 / width TN5 / size 18 - 28
373470 / width TN8 / size 18 - 28



Jumpy

363500 / width TN5 / size 36 - 47



Kilian

374510 / width TN8 / size 23 - 34



Billy

376920 / width TN8 / size 20 - 31



Anden

363600 / width TN5 / size 37 - 47



Monte

363610 / width TN5 / size 21 - 42



Uli

364820 / width TN5 / size 21 - 35
374820 / width TN8 / size 20 - 35



Kerry

364890 / width TN5 / size 25 - 38
374890 / width TN8 / size 20 - 29



Teppi

372580 / width TN8 / size 18 - 35



Hausi

372590 / width TN8 / size 18 - 35



Nils

360030 / width TN5 / size 26 - 33
370030 / width TN8 / size 26 - 34



Kenan

360080 / width TN5 / size 28 - 31
370085 / width TN9 / size 26 - 37



Kito

360750 / width TN5 / size 26 - 35
370750 / width TN8 / size 26 - 29



Fabian

370785 / width TN9 / size 26 - 29



Nico

360730 / width TN5 / size 26 - 38



Kai

362530 / width TN5 / size 18 - 33
372530 / width TN8 / size 18 - 45



Liam

372520 / width TN8 / size 26 - 33
372525 / width TN9 / size 26 - 31



Josephine

362410 / width TN5 / size 20 - 30
372410 / width TN8 / size 20 - 33



Selina

362400 / width TN5 / size 20 - 26
372400 / width TN8 / size 20 - 33

Support shoes

The safe companion – with flexible individual width adjustment

Our aim is to unit a strong team of customers, doctors, physiotherapists, health insurance companies and medical service providers, enabling us to design the best possible therapeutic appliances.

Our support shoes meet all requirements for recognition by health insurance companies as therapeutic appliances. Where each foot requires different care, shoes of differing widths and sizes can be combined without difficulty.

Key shoe quality criteria include:

- Soft materials at the shoe collar in order to avoid pressure points and ensure a high degree of comfort
- Strong heel support while permitting free foot-rolling
- Hook-and-loop fasteners or zips to make shoes easy to put on and take off





Support shoes

Support shoes

Sandals

Page 96–99



Boots

Page 100–109



School shoes

Page 110



Limited quantities

Page 111





Support shoes

Support shoes



Christin

351730 / width TN2



Nicky

351700 / width TN2

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Nicky	351700	TN2	23 - 35	Nubuck	dark blue / black	Hook-and-loop fasteners with roller loop	textile	yes
Christin	351730	TN2	23 - 35	Nubuck	grey / pink	Hook-and-loop fasteners with roller loop	textile	yes



Jack

351020 / width TN4



Tuffy

351090 / width TN4



Sunny

350450 / width TN4

You can use the configurator to design this model for yourself



my design
www.mydesign-schein.de

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Sunny	350450	TN4	23 - 35	Nubuck / Mesh fabric	red / dark blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Jack	351020	TN4	21 - 38	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Tuffy	351090	TN4	30 - 42	Nubuck / Nappa	grey / blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes

Support shoes



You can use the configurator to design this model for yourself



Hero

351085 / width TN2
351080 / width TN4

Domino

351050 / width TN4

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Domino	351050	TN4	27 - 42	Nubuck / Nappa	grey	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Hero	351085 351080	TN2 TN4	23 - 42	Nubuck / Mesh fabric	dark blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes



Support shoes

Support shoes



You can use the configurator to design this model for yourself



Nolly

356350 / width TN2



Yannie

354300 / width TN2



Matti

356390 / width TN2

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Yannie	354300	TN2	23 - 42	Faded denim look	black	Lacing / Zip on the inside	textile	yes
Nolly	356350	TN2	30 - 42	Faded denim look	beige	Lacing / Zip on the inside	textile	yes
Matti	356390	TN2	33 - 45	Nappa	black	Lacing / Zip on the inside	textile	yes



Lani

359410 / width TN2



Jarle

359420 / width TN2



Siva

359430 / width TN2

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Lani	359410	TN2	23 - 39	Nubuck	lilac	Lacing / Zip on the inside	textile	yes
Jarle	359420	TN2	23 - 47	Faded denim look	blue	Lacing / Zip on the inside	textile	yes
Siva	359430	TN2	23 - 35	Faded denim look	pink	Lacing / Zip on the inside	textile	yes

Support shoes



Unique

352680 / width TN2



Glitzi

352660 / width TN2



Momo

358120 / width TN2



my design
www.mydesign-schein.de

You can use the
configurator to design
this model for yourself

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Glitzi	352660	TN2	23 - 32	Nubuck	shiny silver	Hook-and-loop fasteners, top one with roller loop	textile	no
Unique	352680	TN2	27 - 42	Nubuck / Nappa	black / grey	Hook-and-loop fasteners, top one with roller loop	textile	no
Momo	358120	TN2	27 - 35	Nappa	blue	Hook-and-loop fasteners, top one with roller loop	textile	yes



Julian

354170 / width TN2



Ava

354970 / width TN2



River

354980 / width TN2

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Julian	354170	TN2	23 - 42	Nappa	dark blue	Hook-and-loop fasteners, top one with roller loop	textile	yes
Ava	354970	TN2	23 - 35	Nubuck	turquoise	Hook-and-loop fasteners, lower one with no roller loop	textile	no
River	354980	TN2	30 - 47	Nubuck / Nappa	grey / dark grey	Hook-and-loop fasteners, lower one with no roller loop	textile	no

Support shoes



Extra-strong
collar cushioning

Vegas

353420 / width TN2
353425 / width TN4



Phoenix

354480 / width TN2
353480 / width TN4
353485 / width TN6

Stella

353450 / width TN2
353455 / width TN4

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Vegas	353420 353425	TN2 TN4	27 - 42	Nubuck / Mesh fabric	dark blue / grey	Hook-and-loop fasteners with roller loop	textile	yes
Stella	353450 353455	TN2 TN4	27 - 39	Nubuck	purple / grey	Hook-and-loop fasteners with roller loop	textile	yes
Phoenix	354480 353480 353485	TN2 TN4 TN6	23 - 45 23 - 45 23 - 42	Nubuck / Mesh fabric	black / grey	Hook-and-loop fasteners with roller loop	textile	yes



Dani

353460 / width TN2
353465 / width TN4

Support shoes

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Dani	353460 353465	TN2 TN4	23 - 42	Nubuck	blue	Hook-and-loop fasteners with roller loop	textile	yes

Support shoes



Harry

350090 / width TN4
350095 / width TN6



Toby

350770 / width TN4
350775 / width TN6



Matterhorn

353620 / width TN2
353625 / width TN4

Extra-strong
collar cushioning

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Harry	350090 350095	TN4 TN6	23 - 38	Nappa / Mesh fabric	black	Lacing	calf leather lining	no
Toby	350770 350775	TN4 TN6	23 - 38	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Matterhorn	353620 353625	TN2 TN4	30 - 47	Nappa / Mesh fabric	black	Hook-and-loop fasteners, central one with no roller loop	calf leather lining	yes



Support shoes

Support shoes



Teddy fur lining

Nevis

352600 / width TN2



Bobby

352910 / width TN6



Tyler

352900 / width TN6

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Nevis	352600	TN2	23 - 47	Faded-look waxed denim	black	Hook-and-loop fasteners, top one with roller loop	teddy fur	yes
Tyler	352900	TN6	23 - 38	Nappa	black	Lacing	teddy fur	no
Bobby	352910	TN6	23 - 38	Nappa	black	Hook-and-loop fasteners	teddy fur	no



Support shoes

Support shoes

School shoes



Jack

351020 / width TN4 / size 21 - 38
Page 97



Toby

350770 / width TN4 / size 23 - 38
350775 / width TN6 / size 23 - 38
Page 106



Harry

350090 / width TN4 / size 23 - 38
350095 / width TN6 / size 23 - 38
Page 106



Matterhorn

353620 / width TN2 / size 30 - 47
353625 / width TN4 / size 30 - 47
Page 106



Bobby

352910 / width TN6 / size 23 - 38
Page 108



Tyler

352900 / width TN6 / size 23 - 38
Page 108

Support shoes Limited quantities

Limited stocks of these models are available.

Please call us on +49 (0)2191 910 201 to ensure their availability.

Our advice: look for a suitable alternative model from our regular range.



Ben

350215 / width TN6 / size 23 - 42



Pool

351040 / width TN4 / size 31 - 42



Joschi

356310 / width TN2 / size 27 - 39



Happy

356320 / width TN2 / size 27 - 39



Weltmeister

358160 / width TN2 / size 24 - 42



Nele

353495 / width TN6 / size 23 - 33



Mackay

354900 / width TN2 / size 33 - 45

354905 / width TN4 / size 33 - 45



Elena

354930 / width TN2 / size 28 - 33

354935 / width TN4 / size 29 - 35



Fabian

350785 / width TN6 / size 26 - 42



Noah

350675 / width TN6 / size 23 - 34



Carolin

352540 / width TN6 / size 23

Multi-support shoes

Raised upper – plus a sturdy arthrodesis cap

The multi-support shoe has an even more raised upper than our support shoe, plus a sturdy arthrodesis cap. This area can be further reinforced through the use of mouldable thermo-plastic support strips on the inside and outside. These stiffening strips are available in two strengths: strong (white) and medium (skin-coloured). The shoe comes with four white strips. These provide even more lateral and rotational stability in the ankle joint area without hampering the natural rolling of the foot. The integrated rocker bar makes walking easier.

As with our other therapeutic shoes, individual shoe adjustments can be made with these models too.

The applications for the multi-support shoe include the ongoing treatment of ankle joint capsular ligament injuries which originally received conservative or surgical treatment, as well as capsular ligament insufficiency and Achilles tendon injuries. It is also suitable for immobilising the tarsal and mid-foot regions and for the ongoing treatment of fractures and arthrodesis.





Montana

359010 / width ST3



041910-002 Size 36-39
 041910-003 Size 40-44
 041910-004 Size 45-47
 One set contains
 eight white support strips



041911-002 Size 36-39
 041911-003 Size 40-44
 041911-004 Size 45-47
 One set contains
 eight skin-coloured support strips



Nevada

359000 / width ST3



Oregon

359500 / width ST3

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Nevada	359000	ST3	36 - 47	Nappa / Mesh fabric	black / grey	Lacing	textile	no
Montana	359010	ST3	36 - 47	Nappa / Mesh fabric	white	Lacing	textile	no
Oregon	359500	ST3	36 - 47	Nappa / Mesh fabric	black / grey	Hook-and-loop fasteners with roller loop	textile	no

Anti-varus shoes

Anti-varus shoes – for secure correction right from the start

At first, many children find it difficult to adjust to walking in anti-varus shoes. The stiffness of the shoe, which is a must from a therapeutic standpoint, is unaccustomed and makes walking more difficult because the foot is forced into the corrected position. However, experience has shown that children quickly become accustomed to the new feeling.

In cases of slight misalignment, therapy with anti-varus shoes can be supported by the following:

- Regular foot massage
- Suitable exercises
- Physiotherapy

It is important to continuously monitor the treatment progress in consultation with your doctor. The earlier your child receives the right care, the easier it will be to rectify the problem. Where each foot requires different care, shoes of differing widths and sizes can be combined.

Cushioned upper ensures high degree of comfort

High uppers to support the ankle joints

The medial extension of the heel cap beyond the correction point by the big toe's metatarsophalangeal joint ensures consistent correction

Convenient fastening and width adjustment using hook-and-loop straps or laces

Sturdy toe cap to protect toes and prevent wear and tear



Sickle foot and clubfoot orthosis

The sickle foot and clubfoot orthosis for babies is a corrective shoe-shaped orthotic device made of soft leather. To permit variable correction, the sole of the orthosis is divided into two parts. The rear foot is supported by a heel cap, while correction of the forefoot misalignment is done by fitting a tension spring at the other edge of the foot in the Lisfranc joint region. The attachment to the large toe's joint and the cuboid bone area at the outer edge of the foot guides the forefoot into the corrected position. The orthosis is mainly used for infants and toddlers in the crawling age. After that, night splints may be used if further corrective treatment of the forefoot misalignment is necessary.

Pigeon toe orthosis for babies
039101 (right)
039102 (left)

13 sizes from 80 to 150 mm
(not suitable for walking or standing)



Corrective accessory:

The pigeon toe orthosis for babies is delivered with mounted spring at normal tension, plus two springs for medium and maximum tension and a fastening element.

Depending on the prescription,
either dynamic or fixed-position
correction is implemented.



Dynamic correction



Fixed-position correction



Baby sickle foot shoes

The sickle foot shoe for babies is equipped with a fixed degree of anti-varus correction. The shoe's outsole means that it is suitable for standing and walking. In case of one-sided sickle foot, a neutral shoe is available for the unaffected foot.

Eight sizes from 80 to 150 mm

Pigeon toe shoe for babies
with correction
039501 (right)
039502 (left)



Pigeon toe shoe for babies,
neutral position
039506 (right)
039507 (left)



Sickle foot and clubfoot orthosis

The sickle foot and clubfoot orthosis is a passive corrective orthosis made of polyethylene and used to treat childhood sickle foot and previously treated club foot. The orthosis is a night splint to correct the forefoot misalignment by the three-point principle. In case of a slight pointed-foot component, an additional thigh grip can be applied to achieve the desired 90° knee angle, thereby relieving the strain on the calf muscles.

15 sizes
from 7 to 23 cm

Thigh braces
035270 (white)
035271 (blue)
035274 (typhoon)

Lower leg braces
035260 (white)
035261 (blue)
035264 (typhoon)



Anti-varus shoes

Sandals

Page 118–125



Boots

Page 126–141



School shoes

Page 142



Limited quantities

Page 143





Anti-varus shoes

Anti-varus shoes



Amy

041030 / AV
041035 / NS

Lipi

041010 / AV
041015 / NS



Jolina

041050 / AV
041055 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Jolina	041050 041055	AV NS	18 - 29	Patent leather embossed	red	Hook-and-loop fasteners with roller loop	calf leather lining	no
Lipi	041010 041015	AV NS	18 - 42	Nubuck	dark blue / light blue / red	Hook-and-loop fasteners with roller loop	calf leather lining	no
Amy	041030 041035	AV NS	18 - 32	Nubuck	pink / beige	Hook-and-loop fasteners with roller loop	calf leather lining	no



Amelie

341940 / AV
341945 / NS



Kletty

041680 / AV
041685 / NS



Mareike

041620 / AV
041625 / NS



Marina

341910 / AV
341915 / NS

You can use the configurator to design this model for yourself



	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Mareike	041620 041625	AV NS	21 - 32	Nubuck / Mesh fabric	pink	Hook-and-loop fasteners with no roller loop	calf leather lining	yes
Kletty	041680 041685	AV NS	21 - 35	Nubuck / Mesh fabric	dark blue	Hook-and-loop fasteners with no roller loop	calf leather lining	yes
Marina	341910 341915	AV NS	21 - 29	Nubuck	turquoise	Hook-and-loop fasteners with no roller loop	calf leather lining	yes
Amelie	341940 341945	AV NS	21 - 35	Nubuck / Mesh fabric	purple / black	Hook-and-loop fasteners with no roller loop	calf leather lining	yes

Anti-varus shoes



Grey

341880 / AV
341885 / NS



Sofia

341830 / AV
341835 / NS



Anouk

341860 / AV
341865 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Sofia	341830	AV	21 - 29	Patent leather embossed	red	Hook-and-loop fasteners with roller loop	calf leather lining	no
	341835	NS						
Anouk	341860	AV	21 - 29	Nubuck / Nappa	shiny silver / white	Hook-and-loop fasteners with roller loop	calf leather lining	no
	341865	NS						
Grey	341880	AV	21 - 32	Nubuck / Nappa	grey	Hook-and-loop fasteners with roller loop	calf leather lining	no
	341885	NS						



Anti-varus shoes

Anti-varus shoes



Ben

041750 / AV
041755 / NS



Anja

041790 / AV
041795 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Ben	041750 041755	AV NS	18 - 42	Nappa / Nubuck / Mesh fabric	blue / grey / black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Anja	041790 041795	AV NS	18 - 28	Nubuck / Mesh fabric	pink / beige	Hook-and-loop fasteners with roller loop	calf leather lining	no



You can use the
configurator to design
this model for yourself



Mila

341070 / AV
341075 / NS



Emilia

343040 / AV
343045 / NS



Domino

341050 / AV
341055 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Domino	341050 341055	AV NS	21 - 32	Nappa / Nubuck	grey	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Mila	341070 341075	AV NS	18 - 32	Patent leather embossed	red	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Emilia	343040 343045	AV NS	21 - 30	Nubuck	pink	Hook-and-loop fasteners with roller loop	calf leather lining	yes

Anti-varus shoes



Daisy

341140 / AV
341145 / NS



Jack

341020 / AV
341025 / NS



Hero


041080 / AV
041085 / NS

You can use the configurator to design this model for yourself



	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Hero	041080 041085	AV NS	18 - 42	Nubuck / Mesh fabric	dark blue	Hook-and-loop fasteners with roller loop	calf leather lining	yes
Jack	341020 341025	AV NS	21 - 38	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	calf leather lining	no
Daisy	341140 341145	AV NS	21 - 28	Nappa	black	Hook-and-loop fasteners with roller loop	calf leather lining	no

$$E=MC^2$$



Anti-varus shoes



Yellow

349440 / AV
349445 / NS

Siva

349430 / AV
349435 / NS

Lena

349450 / AV
349455 / NS

You can use the configurator to design this model for yourself



Jarle

349420 / AV
349425 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Jarle	349420 349425	AV NS	21 - 42	Faded denim look	blue	Lacing	textile	yes
Siva	349430 349435	AV NS	21 - 29	Denim look	pink	Lacing	textile	yes
Yellow	349440 349445	AV NS	21 - 29	Denim look	yellow	Lacing	textile	yes
Lena	349450 349455	AV NS	21 - 29	Nubuck	turquoise	Lacing	textile	yes



Anti-varus shoes



Yannie

344300 / AV
344305 / NS

You can use the configurator to design this model for yourself



my design
www.mydesign-schein.de



Lina

346360 / AV
346365 / NS



Lara

344330 / AV
344335 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Yannie	344300 344305	AV NS	20 - 42	Faded denim look	black	Lacing	textile	yes
Lara	344330 344335	AV NS	20 - 35	Nappa	pink	Lacing	textile	yes
Lina	346360 346365	AV NS	21 - 35	Nubuck with floral embossing	shiny silver	Lacing	textile	yes



Lynn

040030 / AV
040035 / NS



Louis

340700 / AV
340705 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Lynn	040030 040035	AV NS	21 - 35	Nappa / Nubuck	black	Lacing	calf leather lining	no
Louis	340700 340705	AV NS	23 - 42	Nappa	black	Lacing	calf leather lining	no

Anti-varus shoes



Ava

344970 / AV
344975 / NS



Luana

344910 / AV
344915 / NS



Ole

344920 / AV
344925 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Luana	344910 344915	AV NS	21 - 29	Nubuck	lilac	Hook-and-loop fasteners, lower one with no roller loop	textile	no
Ole	344920 344925	AV NS	21 - 29	Nubuck	royal blue	Hook-and-loop fasteners, lower one with no roller loop	textile	no
Ava	344970 344975	AV NS	21 - 29	Nubuck	turquoise	Hook-and-loop fasteners, lower one with no roller loop	textile	no



Anti-varus shoes



Glitzi

342660 / AV
342665 / NS



Jacko

344180 / AV
344185 / NS



Unique

342680 / AV
342685 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Glitzi	342660 342665	AV NS	21 - 29	Nubuck	shiny silver	Hook-and-loop fasteners, top one with roller loop	textile	no
Unique	342680 342685	AV NS	21 - 39	Nubuck / Nappa	black / grey	Hook-and-loop fasteners, top one with roller loop	textile	no
Jacko	344180 344185	AV NS	23 - 38	Nappa	black	Hook-and-loop fasteners, top one with roller loop	calf leather lining	yes



You can use the configurator to design this model for yourself



Julian

344170 / AV
344175 / NS



Kerry

344890 / AV
344895 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Julian	344170 344175	AV NS	20 - 38	Nappa	dark blue	Hook-and-loop fasteners, top one with roller loop	textile	yes
Momo	348120 348125	AV NS	21 - 32	Nappa	blue	Hook-and-loop fasteners, top one with roller loop	textile	yes
Kerry	344890 344895	AV NS	20 - 38	Nubuck	blue / beige / white	Hook-and-loop fasteners, top one with roller loop	calf leather lining	yes

Anti-varus shoes



Phoenix

042560 / AV
042565 / NS



Vegas

343420 / AV
343425 / NS



Stella

343450 / AV
343455 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Phoenix	042560 042565	AV NS	18 - 42	Nubuck / Mesh fabric	black / grey	Hook-and-loop fasteners with roller loop	textile	yes
Vegas	343420 343425	AV NS	24 - 39	Nubuck / Mesh fabric	dark blue / grey	Hook-and-loop fasteners with roller loop	textile	yes
Stella	343450 343455	AV NS	18 - 32	Nubuck	purple / grey	Hook-and-loop fasteners with roller loop	textile	yes



Extra-strong
collar cushioning

Dani

343460 / AV
343465 / NS



Ida

343430 / AV
343435 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Ida	343430 343435	AV NS	21 - 29	Faded denim look	red	Hook-and-loop fasteners with roller loop	textile	yes
Dani	343460 343465	AV NS	24 - 35	Nubuck	blue	Hook-and-loop fasteners with roller loop	textile	yes

Anti-varus shoes



Extra-strong
collar cushioning


Fabian

040330 / AV
040335 / NS



Toby

340770 / AV
340775 / NS



Harry

340090 / AV
340095 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Fabian	040330 040335	AV NS	18 - 42	Nubuck / Mesh fabric	dark blue / grey	Hook-and-loop fasteners with roller loop	calf leather lining	no
Harry	340090 340095	AV NS	23 - 38	Nappa / Mesh fabric	black	Lacing	calf leather lining	no
Toby	340770 340775	AV NS	23 - 38	Nappa / Mesh fabric	black	Hook-and-loop fasteners with roller loop	calf leather lining	no



Teddy fur
lining



Tyler


342900 / AV
342905 / NS

Bobby

342910 / AV
342915 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Tyler	342900 342905	AV NS	23 - 38	Nappa	black	Lacing	teddy fur	no
Bobby	342910 342915	AV NS	23 - 38	Nappa	black	Hook-and-loop fasteners	teddy fur	no

Anti-varus shoes



Teddy fur
lining

Maxi

042040 / AV
042045 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Maxi	042040 042045	AV NS	18 - 38	Nubuck / Mesh fabric	red / beige / black	Hook-and-loop fasteners with roller loop	teddy fur / SympaTex	yes



Cloud

042030 / AV
042035 / NS



Marty

042050 / AV
042055 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Cloud	042030 042035	AV NS	18 - 42	Nubuck / Mesh fabric	dark blue / black	Hook-and-loop fasteners with roller loop	teddy fur / SympaTex	yes
Marty	042050 042055	AV NS	20 - 42	Nubuck / Mesh fabric	brown / black	Hook-and-loop fasteners with roller loop	teddy fur / SympaTex	yes

Anti-varus shoes

Teddy fur lining



Colima

342630 / AV
342635 / NS



Nevis

342600 / AV
342605 / NS

	Item no.	Width	Size	Material	Colour	Fastener	Inner lining	Reflector
Nevis	342600 342605	AV NS	20 - 42	Faded-look waxed denim	black	Hook-and-loop fasteners, top one with roller loop	teddy fur	yes
Colima	342630 342635	AV NS	21 - 32	Nubuck	purple	Hook-and-loop fasteners with roller loop	teddy fur	yes



Anti-varus shoes

School shoes



Jack

341020 / AV / size 21 - 38
341025 / NS / size 21 - 38
Page 124



Daisy

341140 / AV / size 21 - 38
341145 / NS / size 21 - 38
Page 124



Lynn

040030 / AV / size 21 - 35
040035 / NS / size 21 - 35
Page 129



Louis

340700 / AV / size 23 - 42
340705 / NS / size 23 - 42
Page 129



Jacko

344180 / AV / size 23 - 38
344185 / NS / size 23 - 38
Page 132



Toby

340770 / AV / size 23 - 38
340775 / NS / size 23 - 38
Page 136



Harry

340090 / AV / size 23 - 38
340095 / NS / size 23 - 38
Page 136



Bobby

342910 / AV / size 23 - 38
342915 / NS / size 23 - 38
Page 137



Tyler

342900 / AV / size 23 - 38
342905 / NS / size 23 - 38
Page 137

Anti-varus shoes Limited quantities

Limited stocks of these models are available.

Please call us on +49 (0)2191 910 201 to ensure their availability.

Our advice: look for a suitable alternative model from our regular range.



Jolie

041000 / AV / size 18 - 29
041005 / NS / size 18 - 29



Jamie

041040 / AV / size 18 - 29
041045 / NS / size 18 - 29



Sunny

041090 / AV / size 18 - 35
041095 / NS / size 18 - 35



Linus

341010 / AV / size 18 - 32
341015 / NS / size 18 - 32



Lilly

341030 / AV / size 18 - 32
341035 / NS / size 18 - 30



Sky

041020 / AV / size 18 - 38
041025 / NS / size 18 - 38



Tolouse

341920 / AV / size 21 - 32
341925 / NS / size 22 - 31



Joschi

346310 / AV / size 21 - 32
346315 / NS / size 21 - 32



Happy

346320 / AV / size 21 - 35
346325 / NS / size 21 - 35



Racing

346330 / AV / size 24 - 26
346335 / NS / size 24 - 32



Matti

346390 / AV / size 24 - 36
346395 / NS / size 24 - 36



Annelie

344155 / NS / size 25 - 38



Weltmeister

348160 / AV / size 21 - 30
348165 / NS / size 21 - 32



Teppi

042580 / AV / size 24 - 35
042585 / NS / size 18 - 35



Hausi

042590 / AV / size 24 - 35
042595 / NS / size 18 - 35



Cato

344440 / AV / size 29 - 32
344445 / NS / size 21 - 32



Jona

344420 / AV / size 28 - 32
344425 / NS / size 20 - 32



Kito

040350 / AV / size 32 - 38
040355 / NS / size 20 - 38



Mountain

342610 / AV / size 21 - 32
342615 / NS / size 21 - 32

Orthotic shoes with Boa® fastening system



You can also request our
separate brochure.

We reserve the right to make minor changes to colours and models.
Illustrated shoes are size 28.

Your specialist dealer:



Schein Orthopädie Service KG
Hildegardstr. 5
42897 Remscheid
Germany
Tel. +49 2191 910-201
Fax +49 2191 910-102
remscheid@schein.de
www.mygeneration-schein.de

 carbon neutral
natureOffice.com | DE-221-887193
print production

